




## Accuracy Check Report

Accuracy Check Date	Device ID	Device	User ID	Ambient Temperature
23-Nov-2009 16:09	05254	Pneumotrac	Administrator	23.7°C

Syringe Reference	Syringe Volume	Average Measured Volume	Overall % Difference	Accuracy check result
23390	3 L	3.00 L	0.0%	Pass

**Accuracy Check Results:**



Syringe Stroke	Measured Volume	Duration	% Difference
1	3.00 L	1.32 s	0.0%
2	3.00 L	1.37 s	0.0%
3	3.01 L	3.43 s	0.3%

**Comments:**

Daily Check.

**Signature:** \_\_\_\_\_